



# LENA Q51

Leukemia Fusion Genes DX test



**Make Molecular Detection for  
Blood Disorders Easier, Faster  
and More Comprehensive**



# LENA Q51


Leukemia Fusion Genes DX test

## FDA DISCLAIMER

LENA Q51 ® is offered to professionals only as a "RUO", For research Only and not for Diagnostics purposes.

# Leukemia

Leukemia is an umbrella term for several different cancers of the blood and the blood-forming tissues of the body. All start with problems in the creation of blood cells




# Alercell Provides Comprehensive Solutions for Blood Molecular Detection



Nucleic Acid Extraction System, Lab-Aid 824s  
Total RNA Blood Kit (Lab-Aid)



Nucleic Acid Extraction System, Lab-Aid 896  
Total RNA Blood Kit (Lab-Aid 896)




LENA Q51 Fusion Genes Screening Kit  
Leukemia Fusion Genes Quantification Kit

Sample Extraction

Fluorescent PCR Assay

# Why 52 Genes?



- There are currently more than 200 fusion genes. LENA Q51 actually detects 52 fusion genes. The selection of these 52 fusion genes is to consider an optimal cost performance.
- The 52 fusion genes we selected had a total of more than 200 breakpoints by the time of product development.
- The selection of these 52 fusion genes is based on the fact that the total mutation frequency of these 52 genes exceeds 85%, and the mutation frequency of the remaining more than 100 fusion genes does not exceed 15%.





# LENA Q51



Leukemia Fusion Genes DX test

I

Leukemia Molecular Detection  
Products


# Harm of Leukemia to Human Health Keeps Increasing

- Leukemia, commonly known as blood cancer, is a malignant disease of the blood system
- The incidence rate ranks the first among cancers in children under 15 years old, accounting for about 1/3 of all childhood tumors
- The number of leukemia incidence is increasing year by year, The American Cancer Society's estimates for acute lymphocytic leukemia (ALL) in the United States for 2022 (including both children and adults) are: About 6,660 new cases of ALL (3,740 in males and 2,920 in females)



Source - U.S. Cancer Statistics Working Data by Group CDC. U.S. Cancer Statistics Data Visualizations Tool, based on 2021 submission data (1999- 2019): U.S. Department of Health and Human Services, Centers for Disease Control and Prevention and National Cancer Institute; <https://www.cdc.gov/cancer/dataviz> , released in November 2022.

# More Accurate Classification with Fusion Genes for Precise Diagnosis and Treatment of Leukemia



Easy classification and typing based on FAB (morphology)

**Table 1. WHO classification of myeloid neoplasms and acute leukemia**


**Acute myeloid leukemia (AML) and related neoplasms**

AML with recurrent genetic abnormalities
AML with t(8;21)(q22;q22.1); <i>RUNX1-RUNX1T1</i>
AML with inv(16)(p13.1q22) or t(16;16)(p13.1;q22); <i>CBFB-MYH11</i>
APL with <i>PML-RARA</i>
AML with t(9;11)(p21.3;q23.3); <i>MLLT3-KMT2A</i>
AML with t(6;9)(p23;q34.1); <i>DEK-NUP214</i>
AML with inv(3)(q21.3q26.2) or t(3;3)(q21.3;q26.2); <i>GATA2, MECOM</i>
AML (megakaryoblastic) with t(1;22)(p13.3;q13.3); <i>RBM15-MKL1</i>
<i>Provisional entity: AML with BCR-ABL1</i>
AML with mutated <i>NPM1</i>
AML with biallelic mutations of <i>CEBPA</i>
<i>Provisional entity: AML with mutated RUNX1</i>
AML with myelodysplasia-related changes
Therapy-related myeloid neoplasms

WHO uses **fusion gene** as criteria for leukemia subtype classification since 2016




Adopting next-generation sequencing (NGS) technology for identification of gene fusion in patients with acute leukemia can be a good alternative to conventional tests



(1) Leukemia fusion genes: Acute myeloid leukemia (AML): AML-ETO, PML-RAR $\alpha$ , CBF $\beta$ -MYH11, MLL-AF9 fusion genes. Acute lymphoblastic leukemia (ALL): TEL-AML, E2A-PBX1, BCR-ABL, MLL-AF4 fusion gene. Chronic myeloid leukemia (CML): BCR-ABL fusion gene.


## Fusion Genes are the Most Common Biomarkers in Leukemia



Fusion Genes

- Fusion genes are one of the most common biomarkers in leukemia and are caused by chromosome breakages and recombination.
- More than 200 types of fusion genes have been identified, and about 50% of leukemia patients have fusion genes.
- Different disease subtypes have profiles.

# Fusion Genes Detection Helps Building Precision Leukemia Treatment



## **Fusion Genes Detection Helps Building Precision Leukemia Treatment and is a less expensive path**


- **The average cost of treatment for leukemia patients is \$300,000 to \$900,000**
- **Fusion genes screening charges is less than \$10,000**
- **Quantitative detection charges is even lower four times a year, with an average 10-year survival, total quantitative testing costs less than \$30,000**

**Precise  
medication**

**Timely  
adjustment of the  
treatment  
schedule**

**Relapse warning**

# Provides a Complete Set of Solutions for Leukemia Fusion Genes Detection



# Leukemia Fusion Genes Detection Products



**A complete set of screening and quantification products**

- Leukemia Fusion Genes (Q30) Screening Kit
- Leukemia Fusion Genes (Q51) Screening Kit (Screen the mutations of 30 fusion genes)
- BCR-ABL1 Genotyping Kit (Determine the type: p190, p210 and p230)
- PML-RAR $\alpha$  Genotyping Kit (Determine the type: L, S and V)

- BCR-ABL1 p190 Detection Kit
- BCR-ABL1 p210 Detection Kit
- PML-RAR $\alpha$  L Kit
- PML-RAR $\alpha$  S Kit
- PML-RAR $\alpha$  V Kit
- AML1-ETO Detection Kit
- CBF $\beta$ -MYH11 Kit
- Wt1 Detection Kit

**Clinical Significance**


Drug Prediction

Disease Prognosis Assessment


Treatment Process Monitoring

MRD Detection

# Provides a Complete Set of Solutions for Leukemia Fusion Genes Detection




# Precise Typing on More Fusion Genes to Help Clinical Precision Diagnosis and Treatment

		Competitive product (56)	Clinical significance
PDGFRA Series	<ul style="list-style-type: none"> <li>FIP1L1-PDGFR A</li> <li>ETV6-PDGFR A</li> </ul>	Indistinguishable	<ul style="list-style-type: none"> <li>Seen in about 23% of CELs with a good response to Gleevec</li> <li>Present in AML, tyrosine kinase inhibitors have some effect</li> </ul>
BCR-ABL Series	<ul style="list-style-type: none"> <li>BCR-ABL1 p210</li> <li>BCR-ABL1 p230</li> </ul>	Indistinguishable	<ul style="list-style-type: none"> <li>Seen in about 95% of CML with a good response to Gleevec</li> <li>Rare type, CML or CNL</li> </ul>
AML1 series	<ul style="list-style-type: none"> <li>AML1-MDS1/EV11</li> <li>AML1-MTG16</li> </ul>	Indistinguishable	<ul style="list-style-type: none"> <li><u>1% AML is a poor prognostic factor</u></li> <li><u>Though sensitive to chemotherapy, it has strong toxicity of chemotherapy</u></li> </ul>
MLL series	<ul style="list-style-type: none"> <li>MLL-AF4</li> <li>MLL-AF9</li> <li>MLL-ENL</li> <li>MLL-AF10</li> <li>MLL-SEPT6</li> <li>MLL-ELL</li> <li>MLL-AF17</li> <li>MLL-AF1q</li> <li>MLL-AF1P</li> <li>MLL-AF6</li> <li>MLL-AFX1</li> </ul>	Indistinguishable	<ul style="list-style-type: none"> <li>ALL, poor prognosis, requires intensive therapy</li> <li>AML, poor prognosis, suggests poorer prognosis</li> <li>ALL, mostly neonatal congenital leukemia, poor prognosis</li> <li>AML-M5-type, poor prognosis</li> <li>AML, poor prognosis, slightly better with bone marrow transplantation</li> <li>Adult AML, very poor prognosis, bone marrow transplantation recommended</li> <li>AML</li> <li>AML</li> <li>ALL, AML and MDS, prognosis related to gender and typing</li> <li>AML, very poor prognosis, almost no remission, short survival</li> <li>AML, ALL and CLL in infants and children, poor prognosis</li> </ul>



# Biotech Provides a Complete Set of Solutions for Leukemia Fusion Genes Detection



Philadelphia chromosome (Ph) -like acute lymphoblastic leukemia (ALL), is a high risk subset with a gene expression profile that shares significant overlap with that of Ph -positive (Ph +) ALL

Acute lymphoblastic leukemia (ALL), a malignant clone originating from lymphoid precursor cells, accounts for more than 70% of childhood leukemias.

Currently, treatment schedules based on different subtypes of ALL have achieved good results, but some patients are still poorly treated and prone to relapse. Among them, the newly discovered subtype **Ph-like ALL** is included, accounting for about 10% to 30% of all in adults and 15% of all in children.

**Diagnosis and Treatment Norms for Adult Acute Lymphoblastic Leukemia (2018 Edition)**

Tentative typing: BCR-ABL1 -like ALL

The definition of this type of AL is still difficult. In 2009, two research groups identified a new high-risk subtype of Ph-negative ALL and proposed the concept of Ph-like ALL, which has a similar genetic profile to BCR-ABL 1-positive ALL with IKZF1 or other lymphoid transcriptional regulator deficiency and a similar clinical prognosis, and is therefore called Ph-like ALL (Ph like ALL or BCR/ABL1-like ALL).

WHO lymphoblastic leukemia/lymphoma typing 2016 edition

ALL

B lymphoblastic leukemia/lymphoma

ALL, PTCL-U

ALL with recurrent genetic abnormalities

ALL with t (9; 22) (q34.1; q11.2); BCR-ABL1

ALL with t (v; 11q23.3); KMT2A

ALL with t (12; 21) (p13.2; q22.1) : ETV6-RUNX1

ALL with hyperdiploid karyotype

ALL with hypodiploid karyotype

ALL with t (5; 14) (q31.1; q32.3); IL3-IGH

ALL with t (1; 19) (q23; p13.3); TCF3-PBX1

Tentative typing : BCR-ABL1-like ALL<sup>a</sup>

Tentative typing: B-ALL<sup>a</sup> with internal amplification of chromosome 21





T-lymphoblastic leukemia/lymphoma


Tentative typing: early pre-t-cell lymphoblastic leukemia <sup>a</sup>

Tentative typing: natural killer (NK) cell-lymphocytic leukemia <sup>a</sup>

Note: 1) <sup>a</sup> is a newly added typing. 2) ALL = B lymphoblastic leukemia.

# 4 tubes of Pre-packed Dry Powder Screening Reagent for Selected Ph-like Patients

	Fusion genes	Medication instruction
<b>A tube</b>	 <p><b>ABL1-ETV6</b>, NUP214, RCSD1, RANBP2, SNX2, ZMIZ1, FOXP1 <b>CSF1R-SSBP2</b> <b>IL2RB-MYH9</b></p>	<p>Dasatinib</p> <p><b>JAK2 inhibitors</b> <b>JAK1 and/or inhibitors</b>      JACK3</p>
<b>B tube</b>	 <p><b>JAK2-ATF7IP</b>, BCR, EBF1, ETV6, SSBP2, STRN3, TPR <b>CRLF2-P2RY8</b> <b>TSLP-IQGAP2</b></p>	<p><b>JAK2 inhibitors</b></p> <p><b>JAK2 inhibitors</b> <b>JAK2 inhibitors</b></p>
<b>C tube</b>	 <p><b>JAK2-BCR</b>, PAX5, PPF1BP1, TERF2 <b>NTRK3-ETV6</b> <b>TYK2-MYB</b></p>	<p><b>JAK2inhibitors</b></p> <p>Crizotinib <b>JAK2 inhibitors</b></p>
<b>D tube</b>	 <p><b>ABL2-PAG1</b>, RCSD1, ZC3HAV1 <b>PDGFRB-EBF1</b>, SSBP2, TNIP1, ZEB2, ATF7IP <b>IKZF1</b> (IK6)Independent prognosis</p>	<p>Dasatinib Dasatinib</p>




**Four-tube**, four-color multi-detection system

Pre-packed dry powder form, no solution dispensing required, shipping at room temperature


Detects 32 fusion genes (42 fusion forms) and **one independent prognostic factor**

## Provides a Complete Set of Solutions for Leukemia Fusion Genes Detection





# Quantitative Detection Using an Innovative Single-Tube Dual System, Saving 50% of Manual Operations


Single-tube dual system PK traditional independent system



VS



**New upgrade**  
**Single-tube dual system**



# A Comprehensive and Diverse Range of Quantitative Products to Meet Different Customer Needs



**Operation method:** One-step method  
**Advantages:** detected by adding samples only, easy operation, less contamination, high sensitivity




**Operation method:** Two-step method  
**Advantages:** CDNA samples can be saved and other RNA assays can be performed when detection is carried out

**Two methods are available, with flexible options to meet different clinical needs**

Number	Karyotype	Fusion gene	Number	Karyotype	Fusion gene
1	t(9;11)(p22;q23)	MLL-AF9	31	t(1;11)(q23;p15)	NUP98-PMX1
2	t(15;17)(q24;q21)	PML-RARα	32	t(2;11)(q31;p15)	NUP98-HOXD13
3	t(8;21)(q22;q22)	AML1-ETO	33	t(7;11)(p15;p15)	NUP98-HOXA9
4	t(4;11)(q21;q23)	MLL-AF4	34	t(7;11)(p15;p15)	NUP98-HOXA13
5	t(12;21)(p13;q22)	TEL-AML1	35	t(11;12)(p15;q13)	NUP98-HOXC11
6	t(1;19)(q23;p13)	E2A-PBX1	36	t(7;11)(p15;p15)	NUP98-HOXA11
7	t(11;19)(q23;p13.3)	MLL-ENL	37	der(17)(q21)	STAT5b-RARα
8	t(9;22)(q34;q11)	BCR-ABL1	38	t(11;17)(q13;q21)	NUMA-RARα
9	del(1)(p32)	SIL-TAL1	39	t(4;17)(q12;q21)	FIPIL1-RARα
10	t(10;11)(p12;q23)	MLL-AF10	40	der(17)(q21;q24)	PRKAR1A-RARα
11	inv(16)(p13;q22)	CBFB-MYH11	41	t(10;11)(p13;q21)	CALM-AF10
12	t(3;21)(q26;q22)	AML1-MDS1/EV11	42	t(9;12)(p24;p13)	TEL-JAK2
13	del(4)(q12)	FIPIL1-PDGFRα	43	t(x;11)(q13;q23)	MLL-AFX1
14	del(9)(q34)	SET-CAN	44	t(4;12)(q12;p13)	ETV6-PDGFRα
15	t(17;19)(q22;p13)	E2A-HLF	45	(11q23)	dupMLL
16	t(6;9)(p23;q34)	DEK-CAN	46	t(10;14)(q24;q11)	Hox11
17	t(X;11)(q24;q23)	MLL-SEPT6	47	t(5;14)(q35;q32)	HOX11L2
18	t(16;21)(p11;q22)	TLS-ERG	48	t(3;3)(q21;q26)	EVI-1
19	t(5;12)(q33;p13)	TEL-PDGFRB	49	t(9;22)(q34;q11)	BCR-ABL1 p190
20	t(11;19)(q23;p13.1)	MLL-ELL	50	t(9;22)(q34;q11)	BCR-ABL1 p210
21	t(11;17)(q23;q12-21)	MLL-AF17	51	t(15;17)(q24;q21)	PML-RARα S-type
22	t(5;17)(q35;q21)	NPM-RARα	52	t(15;17)(q24;q21)	PML-RARα V-type
23	t(3;5)(q25;q34)	NPM-MLF1	53	t(15;17)(q24;q21)	PML-RARα L-type
24	t(11;17)(q23;q21)	PLZF-RARα	54	t(2;5)(p23;q35)	NPM1-ALK
25	t(1;11)(q21;q23)	MLL-AF1q			
26	t(1;11)(p32;q23)	MLL-AF1P			
27	t(9;12)(q34;p13)	TEL-ABL1			
28	t(16;21)(q24;q22)	AML1-MTG16			
29	t(3;21)(q26;q22)	AML1-EAP			
30	t(6;11)(q27;q23)	MLL-AF6			

**The most comprehensive fusion gene quantitative detection reagent**

# Easy to Operate and Complete PCR Detection in 2 Hours



Blood total RNA extraction kit and Lab-Aid 824S Nucleic Acid Extraction System, can realize the automatic blood total RNA extraction easily and quickly, subverting all the existing manual extraction methods

# Open Platform for Mainstream Fluorescent PCR Instruments



**ABI 7500**



**Bio-Rad  
CFX96**



**Stratagene  
Mx3000/300  
5P**



**Hongshi SLAN-  
96P/S**

**Applicable for instruments with FAM, HEX/JOE, ROX and Cy5  
detection channels**





# LENA Q51

Leukemia Fusion Genes DX test

## II

Fully Automated Blood RNA  
Extraction System




# Equipped with Own Brand Lab-Aid 824S Automatic Extraction System



- Complete 24 blood RNA extracts in 37 minutes
- No need for pre-treatment, high quality RNA can be obtained in one step of sample addition
- One reagent strip for one sample, flexible for 1 - 24 samples extraction, effectively avoiding reagent waste
- The maximum volume of sample is 2 ml to ensure the total RNA extraction
- Perfectly extract the whole blood samples, leukocyte samples and heparin anticoagulated bone marrow blood samples



# Manual Operation Time Reduced by 90% and Total Time Reduced by 65%



Total extraction time of an imported manual extraction reagent: **150 mins**


Take 24 samples as an example



Total extraction time for Lab-Aid 824s: 52 minutes  
Manual operation : 15 minutes


## High Extraction Efficiency of the Whole Blood Samples, Perfectly Suited for Subsequent Quantitative Assays

Copy Number of Fusion Genes




The copy number of quantitative RNA is comparable between the **200  $\mu$ L** blood sample extracted by Lab-Aid 824s and the **1 mL** blood sample extracted by a competing manual operation, with high extraction efficiency.

# High Quality Extraction of Long-term Preserved Leukocyte Samples



# High Quality Extraction of Heparin-anticoagulated Bone Marrow Blood Samples




■ Heparin anticoagulated bone marrow blood 1  
■ Leukocyte 1  
■ Heparin anticoagulated bone marrow blood 2  
■ Leukocyte 2

--- Standard line

- Lab-Aid 824s extracts 500  $\mu$ L of highly concentrated heparin-anticoagulated bone marrow blood, the RNA concentration, purity and integrity of which can meet the requirements;
- It can remove inhibitors in the blood without affecting subsequent amplification.





- 5 West Mendenhall- Suite 202
- Bozeman MT 59715
- Tel: 888 960 0855
- [Info@Alercell.com](mailto:Info@Alercell.com)